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E-mail: info@ade-hamburg.de

Installation of scale

Your ADE electronic baby hanging scale is a quality product. In order to get the most accurate results please read and respect the instructions of the Operating Manual.

Insert batteries

- Please use two (2) CR2032 (lithium) batteries.
- Open the battery compartment on the back side of the scale by loosening the fixing screws carefully.
- Insert the batteries by respecting the polarity markings (+-).
- Please ensure that the battery compartment cover is being screwed in place before starting a function test.
- When disposing of disused batteries, observe the valid legal and municipal regulations

Display function keys

There are 4 function keys on this scale.

- Key **ON** to start the scale.
- Key **OFF** to switch off the scale.
- Key **SET** to change the weighing mode from kg to lb or to check the room temperature.
- Key **TARE** to reset the scale to zero after weighing procedure.

Operation

Weighing

To achieve a reliable result make sure the scale is placed in a vertical position and no load is hanging at the hook.

3.  Press **ON** key to start the scale. During an automatic zero setting you will see the indication '8888', displayed before to come to '0.00' in kg weighing mode.
4. You may now weigh the good. Make sure to use a weighing bag. The material absolutely has to be a very strong fabric (like canvas) for safety reason. It must be 100% tear proof.

Tare function

Choosing this function you may weigh the baby together with different objects (e.g. nappy, vest etc.) and the weighing bag or – trouser without influencing the weight result.

- Press **ON** and wait for '0.00' indication on the display.
- Weighing bag or trouser including objects (nappy, vest etc.) to be placed on the hook.
- Press **TARE** to reset scale to zero.
- Take all objects away from the scale and press **SET**. The scale will now show only the weight of the weighing bag as negative value. Place the baby incl. objects carefully in the weighing bag or trouser and hang with utmost attention on the hook. Do not swing the baby hanging in the bag.
- The scale will now display the weight of the baby only.
- After the weighing procedure make sure the baby is safely secured in its cradle.
- Press **OFF** to disconnect the scale from power.
- Additionally the scale is equipped with an automatic shut off after 60 seconds without being loaded to prolong the lifetime of the batteries.

Change unit

Should you prefer to weigh in another weighing mode than kg you might change into another mode (lb, oz or indication of room temperature). In this case press **SET** key and choose preferred mode. Make sure the scale has been started before by pressing **ON** and the displays shows '0.00'.

Warning signals and advices

If battery symbol (+-) is indicated on the display, the batteries are too weak to show a proper indication. Replace the batteries.

In case scale will not be used for a long period of time take the batteries out to protect the scale and secure a long lifetime.

If the scale is overloaded (max. capacity 30 kg!) the display shows `over`. Unload the scale.

In case the scale does not show `0.00` after having been started, press TARE to find the zero point.

Cleaning and maintenance

Safety

- Your ADE baby hanging scale is a precision instrument. Handle with care. Do not drop or let heavy object fall onto it.
- Never leave the baby unattended hanging in the weighing trouser.
- Make sure the scale is not getting wet or stored in humid places. Fluids penetrating the scale can severely damage the electronic parts.
- Make sure only discharged batteries respective batteries with isolated poles will be disposed as otherwise short circuits may occur. Never dispose batteries in fire, as they may leak or explode

Cleaning

For cleaning just use a damp cloth or common disinfection products. Do not use any aggressive cleaning fluids. Make sure no water or cleaning fluid penetrates the scale.

Disposal

Electronic devices are no domestic waste.

If the scale is out of use the responsible waste disposal association will be glad to advise you on the measures necessary to dispose of it properly.

In case of faulty results

4. Make sure batteries are inserted correctly.
5. Make sure you are in the weighing mode you selected (kg or lb).
6. Do not load the scale before the display shows `0.00`.

Warranty – Liability

ADE will repair or replace this product within a period of 2 years from date of delivery applying to defects occurred due to poor material or workmanship (presentation of invoice is necessary). All moveable parts – batteries, cables, mains units, rechargeable batteries etc. – are excluded. Warranty does not include natural wear, damages caused by incorrect or negligent handling, overuse or unsuitable or improper storage. The cost of transport will be borne by the customer should the equipment be located anywhere other than the customer's premises. In the event of transport damage, claims under warranty can be accepted subject to using the complete original packaging and the equipment secured the way it was packed and secured originally. All the packaging should therefore be retained.

Any warranty claim will be rejected if the equipment is opened by persons not authorized expressly by ADE to do so.

Customers abroad are kindly requested to contact their local dealer for any warranty claim, customer service and spare parts.

CE Declaration of Conformity

It is hereby certified that this product meets the basic safety requirements laid down in the Directive of the Council applicable to the harmonisation of the legal provisions of the member states with regards to the electromagnetic compatibility (2004/108/EG).

Technical data

- Robust metal housing
- Capacity: 20 kg
- Graduation: 20 g
- Dimensions: 55 x 100 x 22 mm
- LC Display: 15 mm
- Weight: abt 200 g

Supplier

ADE GmbH & Co.
Hammer Steindamm 27 – 29
D – 22089 Hamburg
Germany
☎ +49 (0)40 432776-0
Fax +49 (0)40 432776-10
E-mail: info@ade-hamburg.de